

Superb Nitrile Examination Gloves

Latex Free • Powder Free • Non-Sterile • Disposable • Ambidextrous • Textured

FDA 510(K); ASTM D6319; ASTM D6978; ASTM F1671

Manufacturer: HongRay Group

SUPERB Nitrile gloves combines superb barrier protection, comfort fit, and feel while providing the tensile strength necessary for demanding environments. Our Nitrile gloves are made from 100% Nitrile, a synthetic compound that is completely free of skin-irritating rubber proteins.



Physical Dimensions

Size	Length (mm)	Width (mm)	Thickness (mil)	
			Palm	Finger
S	240±10	85±5	3.5	4
M	240±10	95±5	3.5	4
L	240±10	105±5	3.5	4
XL	240±10	115±5	3.5	4

Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	15	15
Elongation (%)	500	400

Standards and Certifications	ASTM D6319, ASTM D6978, ASTM F1671, Category III, EN 420:2003 + A1:2009, EN 455 1-4, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600, ISO 10993-10:2002, ISO 13485, FDA 510(k) Medical Examination Grade
Freedom from Holes	According to ISO 2859, G-II, multiple sampling plan for normal inspection, AQL 1.5
Powder Residues	For powder free gloves, not more than 2mg per glove
Packing	100 gloves per box; 10 boxes per case; 1000 gloves per case

