

Superb Nitrile Examination Gloves

Latex Free • Powder Free • Non-Sterile • Disposable • Ambidextrous • Textured

ASTM D6319; ASTM D6978; ASTM F1671

Manufacturer

HongRay Group

FDA 510(K) Number

K020493; K182600

Contract Factory (OEM)

Shanxi Hongjin Plastic Technology Co., Ltd

Better Care Plastic Technology Co., Ltd

Factory FDA Registration Number

3013565389

3007842585



Physical Dimensions

Size	Length (mm)	Width (mm)	Thickness (mil)	
			Palm	Finger
S	240±10	85±5	3.5	4
M	240±10	95±5	3.5	4
L	240±10	105±5	3.5	4
XL	240±10	115±5	3.5	4

Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	15	15
Elongation (%)	500	400

Standards and Certifications	ASTM D6319, ASTM D6978, ASTM F1671, Category III, EN 420:2003 + A1:2009, EN 455 1-4, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600, ISO 10993-10:2002, ISO 13485, FDA 510(k) Medical Examination Grade
Freedom from Holes	According to ISO 2859, G-II, multiple sampling plan for normal inspection, AQL 1.5
Powder Residues	For powder free gloves, not more than 2mg per glove
Packing	100 gloves per box; 10 boxes per case; 1000 gloves per case

