

Nitrile Powder Free Examination Gloves

ASTM D6319

FDA 510(K) Number

K153562

Manufacturer

Vietglove Corporation

FDA Registration Number

3012446547



Physical Specifications

Size	Length (mm)(min)	Width (mm)	Thickness (mil)	
			Palm	Finger
XS	240	75±5	3.5±0.8	3.9±0.4
S	240	85±5	3.5±0.8	3.9±0.4
M	240	95±5	3.5±0.8	3.9±0.4
L	240	105±5	3.5±0.8	3.9±0.4
XL	240	115±5	3.5±0.8	3.9±0.4

Item	ASTM D6319		Test Result
	Before Aging (min)	After Accelerated Aging (min)	
Tensile Strength (MPA)	14	14	Pass
Elongation (%)	500	400	Pass

Type	Single Use Non – Sterile Powder Free Nitrile Medical Examination Gloves
Design and Feature	Ambidextrous, straight fingers, finger textured, beaded cuff
Storage Condition	The glove shall maintain their properties when stored in a dry condition. Avoid direct sunlight
Packing	100 gloves per box; 10 boxes per case; 1000 gloves per case

