

Nitrile Gloves

Disposable / Safety / Comfortable / Durable / Powder Free / Textured

ASTM D6319

Material Nitrile (Latex Free)



Physical Dimensions

Size	Length (mm)(min)	Width (mm)	Thickness (mil)(min)	
			Palm	Finger
S	240	85±5	2.8	3.1
M	240	95±5	2.8	3.1
L	240	105±5	2.8	3.1
XL	240	115±5	2.8	3.1

Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	15	14
Elongation (%)	500	400

Freedom from Holes	According to ISO 2859, G-II, multiple sampling plan for normal inspection, AQL 1.5
Powder Residues	For powder free gloves, not more than 2mg per glove

