

# Blue Disposable Nitrile Gloves

**ASTM D6319;**

**Manufacturer**

DN99 International Company Limited

**Factory FDA Registration Number**

3017483450

**Material**

Synthetic Nitrile Latex

**Shelf Life**

3 years from date Manufacture



## Physical Specification (Standard AQL 4.0)

5% deviation is allowed

Size	Length (mm)(min)	Width (mm)	Thickness (mil)	
			Palm	Finger
S	240±2	86±5	2.8	3.1
M	240±2	97±5	2.8	3.1
L	240±2	107±5	2.8	3.1
XL	240±2	115±5	2.8	3.1

Item	ASTM D6319		Test Result
	Before Aging (min)	After Accelerated Aging (min)	
Tensile Strength (MPA)	14	14	Pass
Elongation (%)	500	400	Pass

<b>Type</b>	Non-Sterile Powder Free; Ambidextrous; Finger Tip Textured; Beaded Cuff
<b>Storage</b>	Store in dry and cool Place Temperature not higher than 38°C
<b>Packing</b>	100 gloves per box; 10 boxes per case; 1000 gloves per case

