

Nitrile/Vinyl Blend Examination Gloves Nitrile content min 15%

ASTM D6319; ASTM D6978

Manufacturer

HongRay Group

FDA 510(K) Number

K051660

Contract Factory (OEM)

Grand Work Plastic Products Co., Ltd.

Factory FDA Registration Number

3013565389



Physical Dimensions

Size	Length (mm)(min)	Width (mm)	Thickness (mil)	
			Palm	Finger
S	240	85±5	4.5	5
M	240	95±5	4.5	5
L	240	105±5	4.5	5
XL	240	115±5	4.5	5

Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	16	16
Elongation (%)	450	430

Standards and Certifications	ASTM D6319, ASTM D6978, Category III, EN 420:2003 + A1:2009, EN 455 1-4, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600, ISO 10993-10:2002, ISO 13485, FDA 510(K)
Freedom from Holes	According to ISO 2859, G-II, multiple sampling plan for normal inspection, AQL 2.5
Powder Residues	For powder free gloves, not more than 2mg per glove
Packing	100 gloves per box; 10 boxes per case; 1000 gloves per case

