

Clear Vinyl Examination Gloves

Powder Free • Non-Sterile • Disposable • Ambidextrous • Non-Latex • DOP/DEHP free

FDA 510(K); ASTM D5250

Manufacturer: HongRay Group

Vinyl exam gloves are a cost-effective alternative that are ideal for general patient care



Physical Dimensions

Size	Length (mm)(min)	Width (mm)	Thickness (mil)	
			Palm	Finger
S	230	85±5	3.9	3.9
M	230	95±5	3.9	3.9
L	230	105±5	3.9	3.9
XL	230	115±5	3.9	3.9

Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	11	11
Elongation (%)	350	300

Standards and Certifications	ASTM D5250, EN 420:2003 + A1:2009, FDA21 CFR 177-2600, ISO 13485, US QSR/FDA 510(k) Medical Examination Grade
Freedom from Holes	Inspection Level I - 1.5 AQL
Powder Content	Power-Free
Packing	100 gloves per box; 10 boxes per case; 1000 gloves per case

