

Clear Vinyl Examination Gloves

Powder Free • Non-Sterile • Disposable • Ambidextrous • Non-Latex • DOP/DEHP free

ASTM D5250

Manufacturer

HongRay Group

FDA 510(K) Number

K180381

Contract Factory (OEM)

Shanxi Hongjin Plastic Technology Co., Ltd

Factory FDA Registration Number

3013565389



Physical Dimensions

Size	Length (mm)(min)	Width (mm)	Thickness (mil)	
			Palm	Finger
S	230	85±5	3.9	3.9
M	230	95±5	3.9	3.9
L	230	105±5	3.9	3.9
XL	230	115±5	3.9	3.9

Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	11	11
Elongation (%)	350	300

Standards and Certifications	ASTM D5250, EN 420:2003 + A1:2009, FDA21 CFR 177-2600, ISO 13485, US QSR/FDA 510(k) Medical Examination Grade
Freedom from Holes	Inspection Level I - 1.5 AQL
Powder Content	Power-Free
Packing	100 gloves per box; 10 boxes per case; 1000 gloves per case

