

# Nitrile/Vinyl Blend Gloves

Latex Free • Powder Free • Non-Sterile • Disposable • Ambidextrous

**Brand: Clean Safety**

**Manufacturer: Clean Safety**



## Physical Dimensions

Item#	Size	Glove Weight (Grams)
CS7001	SMALL	5.6
CS7002	MEDIUM	6.0
CS7003	LARGE	6.4
CS7004	X-LARGE	6.8
CS7005	XX-LARGE	Special Order

## Physical Requirements

Packing	100 gloves per box; 10 boxes per case
Tensile Strength	≥ 13 MPA
Elongation	≥300%
AQL	1.5

