

# Nitrile/Vinyl Blend Examination Gloves Nitrile content min 15%

**ASTM D6319; ASTM D6978**

**Brand: Clean Safety**

**Manufacturer: Clean Safety**



## Physical Dimensions

Item#	Size	Average Length (mm)	Width (mm)	Average Thickness (mil)	
				Palm	Finger
CS6001	S	240	85±5	4.5	5
CS6002	M	240	95±5	4.5	5
CS6003	L	240	105±5	4.5	5
CS6004	XL	240	115±5	4.5	5

## Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	16	16
Elongation (%)	450	430

<b>Standards and Certifications</b>	ASTM D6319, ASTM D6978, Category III, EN 420:2003 + A1:2009, EN 455 1-4, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600, ISO 10993-10:2002, ISO 13485, FDA 510(K)
<b>Freedom from Holes</b>	According to ISO 2859, G-II, multiple sampling plan for normal inspection, AQL 2.5
<b>Powder Residues</b>	For powder free gloves, not more than 2mg per glove
<b>Packing</b>	100 gloves per box; 10 boxes per case; 1000 gloves per case

