

Clear Vinyl Examination Gloves

Powder Free • Non-Sterile • Disposable • Ambidextrous • Non-Latex • DOP/DEHP free

FDA 510(K); ASTM D5250

Brand: Clean Safety

Manufacturer: Clean Safety



Vinyl exam gloves are a cost-effective alternative that are ideal for general patient care

Physical Dimensions

Item #	Size	Average Length (mm)	Width (mm)	Average Thickness (mil)	
				Palm	Finger
CS2001	S	240	85±5	4	4
CS2002	M	240	95±5	4	4
CS2003	L	240	105±5	4	4
CS2004	XL	240	115±5	4	4

Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	11	11
Elongation (%)	350	300

Standards and Certifications	ASTM D5250, EN 420:2003 + A1:2009, FDA21 CFR 177-2600, ISO 13485, US QSR/FDA 510(k) Medical Examination Grade
Freedom from Holes	Inspection Level I - 1.5 AQL
Powder Content	Power-Free
Packing	100 gloves per box; 10 boxes per case; 1000 gloves per case

