

Nitrile Disposable Gloves

Latex Free • Powder Free • Non-Sterile • Disposable

Manufacturer

Zhuji Mason Gloves Co., Ltd

FDA Registration Number

3017136893



Physical Dimensions

Size	Length (mm)(min)	Width (mm)	Thickness (mil)	
			Palm	Finger
XS	230	75±3	3.9	4.7
S	230	85±3	3.9	4.7
M	230	95±3	3.9	4.7
L	230	105±3	3.9	4.7
XL	230	115±3	3.9	4.7

Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	14	14
Elongation (%)	500	400

Packing

100 gloves per box; 20 boxes per case; 2000 gloves per case

