

Vinyl Disposable Gloves

Vinyl • Powder free • Latex Free

ASTM D3767; ASTM D5151

Brand: Clean Safety

Manufacturer: Clean Safety

- Ideal for general purpose use and food handling
- Smooth finish for a more natural fit
- Produced without the use of DEHP or DINP



Physical Dimensions

Item#	Size	Average Length (mm)	Width (mm)	Average Thickness (mil)	
				Palm	Finger
CS2101	S	230	85±5	4.5	4.5
CS2102	M	230	95±5	4.5	4.5
CS2103	L	230	105±5	4.5	4.5
CS2104	XL	230	115±5	4.5	4.5

Standards and Certifications	ASTM D3767, ASTM D5151, EN 420, EN 455 FDA21 CFR177-2600, ISO 10993-10:2002, ISO 13485
Freedom from Holes	Inspection Level I - AQL 4.0
Powder Content	Powder-Free
Packing	100 gloves per box; 10 boxes per case; 1000 gloves per case

