

# Clean Safety Superb Nitrile Examination Gloves

Latex Free • Powder Free • Non-Sterile • Disposable • Ambidextrous • Textured

## FDA 510(K); ASTM D6319; ASTM D6978; ASTM F1671

SUPERB Nitrile gloves combines superb barrier protection, comfort fit, and feel while providing the tensile strength necessary for demanding environments. Our Nitrile gloves are made from 100% Nitrile, a synthetic compound that is completely free of skin-irritating rubber proteins.



### Physical Dimensions

| Size | Length (mm) | Width (mm) | Thickness (mil) |        |
|------|-------------|------------|-----------------|--------|
|      |             |            | Palm            | Finger |
| S    | 240±10      | 85±5       | 4               | 4      |
| M    | 240±10      | 95±5       | 4               | 4      |
| L    | 240±10      | 105±5      | 4               | 4      |
| XL   | 240±10      | 115±5      | 4               | 4      |

### Physical Requirements

| Test Result            | Before Aging (min) | After Accelerated Aging (min) |
|------------------------|--------------------|-------------------------------|
| Tensile Strength (MPA) | 15                 | 15                            |
| Elongation (%)         | 500                | 400                           |

|                                     |                                                                                                                                                                                                                   |
|-------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Standards and Certifications</b> | ASTM D6319, ASTM D6978, ASTM F1671, Category III, EN 420:2003 + A1:2009, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600, ISO 10993-10:2002, ISO 13485, FDA 510(k) Medical Examination Grade |
| <b>Freedom from Holes</b>           | According to ISO 2859, G-II, multiple sampling plan for normal inspection, AQL 1.5                                                                                                                                |
| <b>Powder Residues</b>              | For powder free gloves, not more than 2mg per glove                                                                                                                                                               |
| <b>Packing</b>                      | 100 gloves per box; 10 boxes per case; 1000 gloves per case                                                                                                                                                       |

