

TouchFlex Nitrile Examination Glove

Powder Free • Non-Sterile • Disposable • Ambidextrous • Beaded Cuff • Versatile

ASTM D6319;
FDA 510(K) K170136; K110981
Shandong Intco Medical Products Co.,Ltd
Grade Medical
Surface Finger texture
Material NBR



Physical Dimensions

Size	Length (mm) (min)	Width (mm)	Average Thickness (mil)		
			Cuff	Palm	Finger
S	240±10	85±5	2	2.8	3.5
M	240±10	95±5	2	2.8	3.5
L	240±10	105±5	2	2.8	3.5
XL	240±10	115±5	2	2.8	3.5

Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	18	14
Elongation (%)	500	400

Standards and Certifications	ASTM D6319, ISO 9001, ISO 13485, FDA 510(k) Medical Examination Grade
Freedom from Holes	According to ISO 2859, G-I, multiple sampling plan for normal inspection, AQL 1.5
Powder Residues	For powder free gloves, not more than 2mg per glove
Packing	100 gloves per box; 10 boxes per case; 1000 gloves per case

