

# Clean Safety Biodegradable Nitrile Exam Gloves (Black)

Latex Free • Powder Free • Non-Sterile • Disposable • Ambidextrous

- COMPLIED WITH US **FDA** 21 CFR PART 177.2600
- THE PROVEN BIODEGRADABILITY OF A SUBSTANCE HAS BEEN DEMONSTRATED BY UTILIZING THE **ASTM D5511** (REAL-TIME RESULTS)
- CLEARED FOR USE WITH CHEMOTHERAPY DRUGS INCLUDING FENTANYL
- MEDICAL-GRADE NITRILE COMPOSITION

## Physical Dimensions

Size	Length (mm)(min)	Width (mm)	Weight (gram)	Average Thickness (mil)	
				Palm	Finger
M	290	95±10	10.3±0.3	6.7	7.9
L	290	110±10	11.3±0.3	6.7	7.9
XL	290	120±10	12.3±0.3	6.7	7.9
XXL	290	125±10	13.3±0.3	6.7	7.9

## Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	14	14
Elongation (%)	500	400

Freedom from Holes	According to ISO 2859, G-II, multiple sampling plan for normal inspection, AQL 1.5
Powder Residues	For powder free gloves, not more than 2mg per glove

