

Superb Black Nitrile Examination Gloves

Latex Free • Powder Free • Non-Sterile • Disposable • Ambidextrous

CleanSafety Superb Black Nitrile gloves Provide Superb protection and performance combined with superb comfort and feel that offers effective barrier protection and tensile strength.

Behind each nitrile examination glove is the assurance of our Clean Safety team. We conduct meticulous quality control measures while aligning with **ASTM**, **ISO**, and **FDA** standards to ensure the performance and comfort of our products. In turn, we proudly offer a premium alternative to latex and vinyl gloves alike.



Physical Dimensions

Size	Product Code	Length (mm/inches)	Thickness (mil)	
			Palm	Finger
S	CS4901	230 / 9.1	2.8	3
M	CS4902	230 / 9.1	2.8	3
L	CS4903	230 / 9.1	2.8	3
XL	CS4904	230 / 9.1	2.8	3

Features	Textured fingertips for added grip Composed of 100% nitrile Medical-grade nitrile composition
Benefits	Lightweight and flexible for added comfort Excellent sensitivity and strength 100 gloves per dispenser for best-in-class savings Versatile applications across industries
General Specs	Nitrile, Powder Free, Textured, 9" length, Black, Healthcare/Industrial/Food Service use

