

Clean Safety Biodegradable Nitrile Exam Gloves (Blue)

Latex Free • Powder Free • Non-Sterile • Disposable • Ambidextrous

- COMPLIED WITH US **FDA** 21 CFR PART 177.2600
- THE PROVEN BIODEGRADABILITY OF A SUBSTANCE HAS BEEN DEMONSTRATED BY UTILIZING THE **ASTM D5511** (REAL-TIME RESULTS)
- CLEARED FOR USE WITH CHEMOTHERAPY DRUGS INCLUDING FENTANYL
- MEDICAL-GRADE NITRILE COMPOSITION

Physical Dimensions

Size	Length (mm)(min)	Width (mm)	Weight (gram)	Average Thickness (mil)	
				Palm	Finger
S	230	80±10	3.1±0.1	3.5±0.5	4.0±0.5
M	230	95±10	3.4±0.1	3.5±0.5	4.0±0.5
L	230	110±10	3.7±0.1	3.5±0.5	4.0±0.5
XL	230	120±10	3.9±0.1	3.5±0.5	4.0±0.5

Physical Requirements

Test Result	Before Aging (min)	After Accelerated Aging (min)
Tensile Strength (MPA)	14	14
Elongation (%)	500	400

Freedom from Holes	According to ISO 2859, G-II, multiple sampling plan for normal inspection, AQL 1.5
Powder Residues	For powder free gloves, not more than 2mg per glove

